Case 1 Neutrophilia

A woman in her 90s undergoing chemotherapy for diffuse large cell lymphoma.

Blood smear (May-Giemsa staining)







PB (×1,000)



Visual differential counts

Blast 0.0 Promyelo 0.0/ 0.0/ Myelo 2.0/ Meta 43.0 Band Seg 53.0 0.0 Eosino Baso 0.0 Mono 1.0 1.0 Lympho At-Ly 0.0 NRBC Other

Celltac Data

Numerical results

WBC	24.33	Н	10³/µl
RBC	4.40		$10^{6}/\mu$ l
HGB	12.42		g/dL
НСТ	37.2		%
MCV	84.5		fL
MCH	28.2		pg
MCHC	33.4		g/dL
RDW-CV	13.4		%
RDW-SD	45.3		fL
PLT	274.8		10³/µl
PCT	0.21		%
MPV	7.8		fL
PDW	16.4		%
P-LCR	32.9		%
P-LCC	90.4		10³/µl
NE	23.19	*	10³/µl
LY	0.62	*	10³/µl
MO	0.20	*	10³/µl
EO	0.26	*	10³/µl
BA	0.06	*	10³/µl
NE%	95.32	*	%
LY%	2.54	*	%
MO%	0.84	*	%
EO%	1.05	*	%
BA%	0.25	*	%



WBC Count

Flags

Morphological Flags	Numerical Flags
Blast Immature Granulocyte	Leukocytosis Neutrophilia
Left Shift Atvoical Ly	Lymphopenia
Ly-Mo Interference	

Explanation of scattergram / histogram

Many plots appear in the neutrophil area (O) on the MAIN scattergram, suggesting neutrophilia. A "neutrophilia" flag is displayed to indicate this. Additionally, in the NE-EO scattergram, the neutrophil plot shows a distribution that extends to the top (O), and the "immature Gr" and "left shift" flags are displayed, suggesting the appearance of immature granulocytes and an increase in band neutrophils.

Explanation of case

A complete blood count revealed an elevated number of white blood cells. A visual white blood cell differential count showed elevated neutrophil levels with left shift, with 2.0 % metamyelocytes, 43.0 % band neutrophils, and 53.0 % segmented neutrophils. This is considered to have been due to the administration of G-CSF during chemotherapy. Smear samples indicated echinocytes. Could this have been due to the infusion? Many band neutrophils with some vacuoles were observed, but no toxic granules or Döhle bodies were identified.

Data provision and supervision, Department of Clinical Laboratory, The University of Tokyo Hospital Masahiro Jona



Scattergrams



Histograms



Case 2 Eosinophilia

A woman in her 50s visited the hospital for atopic dermatitis.

Blood smear (May-Giemsa staining)







PB (×1,000)



Visual differential counts

0.0 Blast Promyelo 0.0 0.0/ Myelo 0.0/ Meta Band 0.0 Seg 24.0 43.0 Eosino Baso 2.0 Mono 4.0 27.0 Lympho At-Ly 0.0 NRBC 0.0 Other

Flags

Ne-Eo Interference

Celltac Data

10.65

4.19

13.03

37.7

90.0

31.1

34.6

12.2

43.9

333.5

0.32

9.5

17.1

51.5

171.8

2.49

2.56

0.40

5.12

0.08

23.37

24.05

3.74

48.08

0.76

0.80

6.4

93.6

5.1

1.3

0.0335

н

Н

Numerical results

WBC

RBC

HGB

HCT

MCV

MCH

MCHC

RDW-CV

RDW-SD

PLT

PCT

MPV

PDW

P-LCR

P-LCC

NE

LY

MO

ΕO

BA

NE%

LY%

MO%

EO%

BA%

RET

IRF

LFR

MFR

HFR

RET%





Explanation of case

A complete blood count revealed an elevated number of white blood cells. A visual white blood cell differential count indicated 43.0 % eosinophils. No other abnormal findings were observed. Gene mutations related to eosinophilia, such as BCR::ABL, PDGFRA, PDGFRB, and FGFR1, were negative. Hence, eosinophilia due to atopic dermatitis, were considered. Smear samples showed an increase in mature eosinophils, which were approximately 16 µm in size and twice the size of the red blood cells.

Explanation of scattergram / histogram

Many eosinophil plots can be observed in the MAIN scattergram and NE-EO scattergram, suggesting eosinophilia (O). A "eosinophilia" flag is displayed to indicate this.

10³/µL

10⁶/µL

g/dL

%

fl

pg

%

fL

%

fL

%

%

10³/µL

10³/µL

10³/µL

10³/µL

10³/µL

10³/µL

%

%

%

%

%

%

%

%

%

%

10⁶/µL

10³/µL

g/dL



WBC Count

FL525/FSC





Scattergrams





Histograms



RET Scattergrams



Data provision and supervision, Department of Clinical Laboratory, The University of Tokyo Hospital Masahiro Jona

Case 3 Basophilia

A woman in her 50s undergoing follow-up observation for polycythemia vera (JAK2V617F positive) without treatment.

Blood smear (May-Giemsa staining)





PB (× 400)



PB (×1,000)



PB (×1,000)

Explanation of case

A complete blood count revealed an elevated number of white blood cells, red blood cells, and platelets. A visual white blood cell differential count showed an increase in basophils to 3.0 %. This was considered to be due to the disease. Smear samples indicated an increase in mature basophils. Basophil granules are water-soluble, so some basophils have fewer granules. Observations of scattered large platelets were made.

Visual differential counts

Blast 0.0 Promyelo 0.0/ 0.0/ Myelo Meta 0.0/ 1.5/ Band Seg /81.5/ 5.5 Eosino Baso 3.0 Mono 1.5 7.0 Lympho At-Ly 0.0 NRBC 0.0 Other

Celltac Data

Numerical results

VBC	34.07	Н	10³/µl
RBC	7.41	Н	10 ⁶ /µl
IGB	15.03		g/dL
ICT	49.0		%
ICV	66.1	L	fL
ICH	20.3	L	pg
ICHC	30.7	L	g/dL
RDW-CV	17.6	Н	%
RDW-SD	46.5		fL
νLT	849.9	Н	10³/µ
СТ	0.71	Н	%
1PV	8.4		fL
WD	20.1	Н	%
-LCR	38.8		%
P-LCC	329.8	Н	10³/µl
١E	25.47	*	10³/µl
Y	2.75	*	10³/µl
10	1.50	*	10³/µl
0	3.11	*	10³/µl
BA	1.24	*	10³/µ
IE%	74.73	*	%
Y%	8.08	*	%
10%	4.41	*	%
0%	9.14	*	%
BA%	3.64	*	%
RET	0.1252		10º/µl
RET%	1.69		%
RF	14.1		%
FR	85.9	L	%
1FR	9.2		%
IFR	4.9	Н	%

Size





RNF FL525/FSC

Flags

Morphological Flags Numerical Flags Blast Immature Granulocyte Left Shift Ne-Eo Interference

Leukocytosis Neutrophilia Monocytosis Eosinophilia Basophilia Erythrocytosis Microcvtosis Thrombocytosis

Explanation of scattergram / histogram

Many plots observed in the basophil region (O) on the MO-BA scattergram, suggesting basophilia. A "basophilia" flag is displayed to indicate this.



Scattergrams

Histograms



RET Scattergrams



Data provision and supervision, Department of Clinical Laboratory, The University of Tokyo Hospital Masahiro Jona

Case 4

AML, myelodysplasia related(AML-MR)

A woman in her 70s visited the hospital with essential thrombocythemia (CALR mutation positive). HU (hydroxyurea) was prescribed, and 10 years later, blast cells were observed. A diagnosis of MDS-IB1 was made, together with follow-up observations.

Blood smear (May-Giemsa staining)





Visual differential counts

78.0

0.0

0.0

1.0

1.0 11.0

0.0

0.0

1.0

8.0

0.0

0.0

0.0

Blast

Myelo

Meta Band

Seg

Eosino

Baso

Mono

At-Ly

NRBC

Other

Lympho

Promyelo





PB (×1,000)





Explanation of case

A complete blood count revealed an elevated number of white blood cells, anemia, and thrombocytopenia. A visual white blood cell differential count revealed 78.0% blast cells. Dysplasia of degranulated neutrophils was also observed in the neutrophils. Due to the history of chemotherapy treatment, it was thought that the case involved a progression to treatment-related myeloid neoplasm (AML-MR). Smear samples showed a cell size of 15–20 μ m, basophilic cytoplasm, delicate nuclear reticulum, nucleolus and blast cells with cleft nucleus. Degranulated neutrophils and hyper segmented neutrophils with degranulation were also observed.

Celltac Data

Numerical results

WBC	11.60	Н	10³/µ
RBC	2.89	L	$10^{6}/\mu$
HGB	8.65	L	g/dL
HCT	25.3	L	%
MCV	87.5		fL
MCH	29.9		pg
MCHC	34.2		g/dL
RDW-CV	12.8		%
RDW-SD	44.8		fL
PLT	20.7	*	10³/µ
PCT	0.02	L	%
MPV	9.2		fL
PDW	20.4	Н	%
P-LCR	50.4		%
P-LCC	10.4	L	10³/µ
NE	0.80	*	10³/µ
LY	4.89	*	10³/µ
MO	5.33	*	10³/µ
EO	0.00	*	10³/µ
BA	0.58	*	10³/µ
NE%	6.91	*	%
LY%	42.13	*	%
MO%	45.94	*	%
EO%	0.00	*	%
BA%	5.02	*	%
RET	0.0020	L	10 ⁶ /µ
RET%	0.07	L	%
IRF	14.3		%
LFR	85.7	L	%
MFR	14.3		%
HFR	0.0		%



Histograms

ľ	TDC
С	ount
	2
E.	

Flags

Morphological Flags	Numerical Flag
Blast	Neutropenia
Left Shift	Lymphocytosis
Atypical Ly	Monocytosis
Ly-Mo Interference	Basophilia
	Anemia
	Thrombocytopeni



Explanation of scattergram / histogram

Monocyte plots on the MAIN scattergram and MO-BA scattergram show an abnormal distribution extending to the top (O) and distributed even in the blast cell flag detection area, suggesting the appearance of blast cells. A "blast" flag is displayed to indicate this.

The reticulocyte plot is almost invisible in the RET scattergram, indicating a decrease in the number of reticulocytes.

Data provision and supervision, Department of Clinical Laboratory, The University of Tokyo Hospital Masahiro Jona

L525/FSC







Granularity

RET Scattergrams

