

**MEK-3DN** HEMATOLOGY CONTROL 2.0mL

ASSAY SHEET

LOT	B242N
	2024-05-05
OPEN	14DAYS

**MEK-1301,1302,1305 Barcode Entry**

The barcodes may be read in any order. If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420	MEK -6500 -6510	MEK -6450	MEK -6550	MEK -1301 -1302 -1305
WBC	10 <sup>3</sup> /μL	±0.8	7.9	7.9	7.9	7.9	7.75
RBC	10 <sup>6</sup> /μL	±0.20	4.66	4.66	4.66	4.66	4.66
HGB	g/dL	±0.5	13.9	14.1	13.9	14.1	13.95
HCT	%	±3.0	40.0	40.0	40.0	40.0	45.9
MCV	fL	±7.0	85.9	85.9	85.9	85.9	98.5
MCH	pg	±2.4	29.8	30.3	29.8	30.3	29.9
MCHC	g/dL	±3.0	34.8	35.3	34.8	35.3	30.4
PLT	10 <sup>3</sup> /μL	±50	278	278	278	278	281.8
MPV	fL	±2.0	5.8	8.2	5.8	8.2	8.4
RDW-CV *	%	±3.0	16.3	16.3	16.3	16.3	15.6
LY%	%	±15.0	33.8	33.8	33.8	33.8	34.64
MO%	%	±8.0	5.4	5.4	5.4	5.4	2.69
GR%	%	±15.0	60.8	60.8	60.8	60.8	62.67
LY	10 <sup>3</sup> /μL	±1.5	2.7	2.7	2.7	2.7	2.68
MO	10 <sup>3</sup> /μL	±0.8	0.4	0.4	0.4	0.4	0.21
GR	10 <sup>3</sup> /μL	±1.5	4.8	4.8	4.8	4.8	4.86

<b>MEK-1305</b>
QC ESR
10±5 mm

<b>Manual Registration</b>	
Lot	242
Level	Normal
Control	MEK-3D

\* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.



<https://www.nihonkohden.com/spi/index.html>



**NIHON KOHDEN CORPORATION**  
1-31-4 NISHIOCHIAI, SHINJUKU-KU  
TOKYO 161-8560 JAPAN



0694-9016850

**MEK-3DH HEMATOLOGY CONTROL 2.0mL**

**ASSAY SHEET**

<b>LOT</b>	B242H
	2024-05-05
<b>OPEN</b>	14DAYS

**MEK-1301,1302,1305 Barcode Entry**

The barcodes may be read in any order. If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.



		RANGE	MEK -6400 -6410 -6420	MEK -6500 -6510	MEK -6450	MEK -6550	MEK -1301 -1302 -1305
<b>WBC</b>	10 <sup>3</sup> /μL	±2.2	20.9	20.9	20.9	20.9	20.47
<b>RBC</b>	10 <sup>6</sup> /μL	±0.25	5.87	5.87	5.87	5.87	5.90
<b>HGB</b>	g/dL	±0.7	19.3	19.3	19.3	19.3	19.42
<b>HCT</b>	%	±4.0	54.1	54.1	54.1	54.1	62.6
<b>MCV</b>	fL	±7.0	92.2	92.2	92.2	92.2	106.1
<b>MCH</b>	pg	±2.8	32.9	32.9	32.9	32.9	32.9
<b>MCHC</b>	g/dL	±3.0	35.7	35.7	35.7	35.7	31.0
<b>PLT</b>	10 <sup>3</sup> /μL	±70	592	592	592	592	583.7
<b>MPV</b>	fL	±3.0	5.8	8.1	5.8	8.1	8.5
<b>RDW-CV *</b>	%	±3.0	16.2	16.2	16.2	16.2	17.1
<b>LY%</b>	%	±15.0	19.6	19.6	19.6	19.6	14.75
<b>MO%</b>	%	±8.0	3.3	3.3	3.3	3.3	1.65
<b>GR%</b>	%	±15.0	77.1	77.1	77.1	77.1	83.60
<b>LY</b>	10 <sup>3</sup> /μL	±3.0	4.1	4.1	4.1	4.1	3.02
<b>MO</b>	10 <sup>3</sup> /μL	±1.5	0.7	0.7	0.7	0.7	0.34
<b>GR</b>	10 <sup>3</sup> /μL	±3.0	16.1	16.1	16.1	16.1	17.11

Manual Registration	
Lot	242
Level	High
Control	MEK-3D

\* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.



<https://www.nihonkohden.com/spi/index.html>




**NIHON KOHDEN CORPORATION**  
1-31-4 NISHIOCHIAI, SHINJUKU-KU  
TOKYO 161-8560 JAPAN



0694-9016860

**MEK-3DL**
**HEMATOLOGY CONTROL**
**2.0mL**
**ASSAY SHEET**

<b>LOT</b>	B242L
	2024-05-05
<b>OPEN</b>	14DAYS

**MEK-1301,1302,1305 Barcode Entry**

The barcodes may be read in any order.  
If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420	MEK -6500 -6510	MEK -6450	MEK -6550	MEK -1301 -1302 -1305
<b>WBC</b>	10 <sup>3</sup> /μL	±0.6	2.0	2.0	2.0	2.0	2.02
<b>RBC</b>	10 <sup>6</sup> /μL	±0.15	2.41	2.41	2.41	2.41	2.36
<b>HGB</b>	g/dL	±0.4	6.0	6.1	6.0	6.1	6.06
<b>HCT</b>	%	±2.0	17.7	17.7	17.7	17.7	19.6
<b>MCV</b>	fL	±9.0	73.5	73.5	73.5	73.5	83.0
<b>MCH</b>	pg	±2.4	24.9	25.3	24.9	25.3	25.7
<b>MCHC</b>	g/dL	±3.0	33.9	34.5	33.9	34.5	30.9
<b>PLT</b>	10 <sup>3</sup> /μL	±30	67	67	67	67	75.0
<b>MPV</b>	fL	±3.0	5.8	8.2	5.8	8.2	8.1
<b>RDW-CV *</b>	%	±3.0	18.2	18.2	18.2	18.2	16.7
<b>LY%</b>	%	±10.0	68.6	68.6	68.6	68.6	68.54
<b>MO%</b>	%	±5.0	6.9	6.9	6.9	6.9	3.94
<b>GR%</b>	%	±10.0	24.5	24.5	24.5	24.5	27.52
<b>LY</b>	10 <sup>3</sup> /μL	±0.5	1.4	1.4	1.4	1.4	1.38
<b>MO</b>	10 <sup>3</sup> /μL	±0.2	0.1	0.1	0.1	0.1	0.08
<b>GR</b>	10 <sup>3</sup> /μL	±0.5	0.5	0.5	0.5	0.5	0.56

**MEK-1305**
**QC ESR**
**48±17 mm**
**Manual Registration**

Lot	242
Level	Low
Control	MEK-3D

\* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.


<https://www.nihonkohden.com/spi/index.html>

**NIHON KOHDEN CORPORATION**

1-31-4 NISHIOCHIAI, SHINJUKU-KU

TOKYO 161-8560 JAPAN



MEK-3DL



B242L



2024-05-05



0694-9016840