


MK-RE2 HEMATOLOGY CONTROL
for reticulocyte 3.0mL

ASSAY SHEET

LOT	MKRE245-2
	2024-07-05
OPEN	14 DAYS

		MEK-9200	RANGE
RET%	%	4.15	±2.85
RET	10 ⁶ /μL	0.1488	±0.1022

2D Barcode Entry



• If RBC value is under 3.00 [$\times 10^6/\mu\text{L}$], please refrain from using.
It may have deteriorated or hemolyzed, which affects the measurement accuracy.



<https://www.nihonkohden.com/spi/index.html>

Manual Registration						
Lot	245					
Level	2					
Expiration	7	M	5	D	2024	Y
Control	MKRE					




NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN

MK-RE3
HEMATOLOGY CONTROL
 for reticulocyte

3.0mL

ASSAY SHEET

LOT	MKRE245-3
	2024-07-05
OPEN	14 DAYS

		MEK-9200	RANGE
RET%	%	8.12	±3.88
RET	10⁶/μL	0.2926	±0.1398

2D Barcode Entry


• If RBC value is under 3.00 [$\times 10^6/\mu\text{L}$], please refrain from using.
 It may have deteriorated or hemolyzed, which affects the measurement accuracy.


<https://www.nihonkohden.com/spi/index.html>
Manual Registration

Lot	245					
Level	3					
Expiration	7	M	5	D	2024	Y
Control	MKRE					



 0123



NIHON KOHDEN CORPORATION
 1-31-4 NISHIOCHIAI, SHINJUKU-KU
 TOKYO 161-8560 JAPAN

MK-RE1 HEMATOLOGY CONTROL
for reticulocyte **3.0mL**

ASSAY SHEET

LOT	MKRE245-1
	2024-07-05
OPEN	14 DAYS

		MEK-9200	RANGE
RET%	%	1.56	±1.44
RET	10 ⁶ /μL	0.0556	±0.0512

2D Barcode Entry



• If RBC value is under 3.00 [$\times 10^6/\mu\text{L}$], please refrain from using.
It may have deteriorated or hemolyzed, which affects the measurement accuracy.



<https://www.nihonkohden.com/spi/index.html>

Manual Registration						
Lot	245					
Level	1					
Expiration	7	M	5	D	2024	Y
Control	MKRE					



NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN