


**MEK-CAL** HEMATOLOGY CALIBRATOR 2.0mL

**ASSAY SHEET**

<b>LOT</b>	PLUS252
	2025-03-05
<b>OPEN</b>	7DAYS

MEK-9100,9200  
2D Barcode Entry



MEK-1301,1302,1303,1305  
Barcode Entry



		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305
WBC	10 <sup>3</sup> /μL	9.3	9.3	9.6	9.5	9.3	9.69	9.30
RBC	10 <sup>6</sup> /μL	4.70	4.70	4.68	4.70	4.67	4.68	4.70
HGB	g/dL	13.7	13.8	13.6	13.5	13.4	13.75	13.66
HCT	%	40.0	40.0	41.9	40.7	39.6	46.1	46.4
MCV	fL	85.2	85.2	89.5	86.5	84.7	98.4	98.7
MCH	pg	29.1	29.4	29.1	28.7	28.7	29.4	29.1
MCHC	g/dL	34.3	34.5	32.5	33.2	33.8	29.8	29.4
PLT	10 <sup>3</sup> /μL	242	242	234	237	236	253.2	263.0
MPV	fL	6.4	8.6	6.4	7.4	6.4	9.0	9.0
RDW-CV *	%	15.5	15.5	14.5	15.6	14.7	15.6	15.5

MEK	-7222
	-7300
	-8222
	-9100
	-9200
FS	168
FL	87
SD	105
FS THR.	100

\* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,7222,8222.



<https://www.nihonkohden.com/spi/index.html>

**NIHON KOHDEN CORPORATION**  
1-31-4 NISHIOCHIAI, SHINJUKU-KU  
TOKYO 161-8560 JAPAN

