


**MK-RE2**

**HEMATOLOGY CONTROL**  
for reticulocyte

3.0mL

**ASSAY SHEET**

LOT	MKRE253-2
	2025-05-05
OPEN	14 DAYS

		MEK-9200	RANGE
RET%	%	4.46	±2.54
RET	10 <sup>6</sup> /μL	0.1553	±0.0885

**2D Barcode Entry**



- If RBC value is under 3.00 [ $\times 10^6/\mu\text{L}$ ], please refrain from using.  
It may have deteriorated or hemolyzed, which affects the measurement accuracy.



<https://www.nihonkohden.com/spi/index.html>

**Manual Registration**

Lot	253					
Level	2					
Expiration	5	M	5	D	2025	Y
Control	MKRE					




**NIHON KOHDEN CORPORATION**  
1-31-4 NISHIOCHIAI, SHINJUKU-KU  
TOKYO 161-8560 JAPAN

**MK-RE3**

**HEMATOLOGY CONTROL**  
for reticulocyte

3.0mL

**ASSAY SHEET**

LOT	MKRE253-3
	2025-05-05
OPEN	14 DAYS

		MEK-9200	RANGE
RET%	%	8.70	±3.30
RET	10 <sup>6</sup> /μL	0.3087	±0.1171

**2D Barcode Entry**



- If RBC value is under 3.00 [ $\times 10^6/\mu\text{L}$ ], please refrain from using.  
It may have deteriorated or hemolyzed, which affects the measurement accuracy.



<https://www.nihonkohden.com/spi/index.html>

**Manual Registration**

Lot	253					
Level	3					
Expiration	5	M	5	D	2025	Y
Control	MKRE					




**NIHON KOHDEN CORPORATION**  
1-31-4 NISHIOCHIAI, SHINJUKU-KU  
TOKYO 161-8560 JAPAN

**MK-RE1**
**HEMATOLOGY CONTROL**  
 for reticulocyte

3.0mL

**ASSAY SHEET**

LOT	MKRE253-1
	2025-05-05
OPEN	14 DAYS

		MEK-9200	RANGE
RET%	%	1.56	±1.44
RET	10 <sup>6</sup> /μL	0.0545	±0.0504

**2D Barcode Entry**


- If RBC value is under 3.00 [ $\times 10^6/\mu\text{L}$ ], please refrain from using.  
It may have deteriorated or hemolyzed, which affects the measurement accuracy.


<https://www.nihonkohden.com/spi/index.html>
**Manual Registration**

Lot	253					
Level	1					
Expiration	5	M	5	D	2025	Y
Control	MKRE					


**NIHON KOHDEN CORPORATION**  
 1-31-4 NISHIOCHIAI, SHINJUKU-KU  
 TOKYO 161-8560 JAPAN