

MEK-3DN

HEMATOLOGY CONTROL

2.0mL

ASSAY SHEET

LOT	B255N
	2025-08-05
OPEN	14DAYS

MEK-1301,1302,1305 Barcode Entry

The barcodes may be read in any order.

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
WBC	10 ³ /μL	±0.8	8.0	8.0	7.82
RBC	10 ⁶ /μL	±0.20	4.89	4.89	4.91
HGB	g/dL	±0.5	14.1	14.2	14.11
HCT	%	±3.0	41.6	41.6	47.5
MCV	fL	±7.0	85.0	85.0	96.7
MCH	pg	±2.4	28.8	29.0	28.7
MCHC	g/dL	±3.0	33.9	34.1	29.7
PLT	10 ³ /μL	±50	251	251	274.6
MPV	fL	±2.0	5.3	8.0	8.2
RDW-CV *	%	±3.0	16.7	16.7	17.1
LY%	%	±10.0	33.0	33.0	30.49
MO%	%	±5.0	3.4	3.4	2.55
GR%	%	±10.0	63.6	63.6	66.96
LY	10 ³ /μL	±1.0	2.6	2.6	2.38
MO	10 ³ /μL	±0.4	0.3	0.3	0.20
GR	10 ³ /μL	±1.0	5.1	5.1	5.24

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.

MEK-1305

QC ESR

8±5 mm

Manual Registration

Lot	255
Level	Normal
Control	MEK-3D


<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU

TOKYO 161-8560 JAPAN



MEK-3DN



B255N



2025-08-05



0694-901685T

MEK-3DH

HEMATOLOGY CONTROL

2.0mL

ASSAY SHEET

LOT	B255H
	2025-08-05
OPEN	14DAYS

MEK-1301,1302,1305 Barcode Entry

The barcodes may be read in any order.

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
WBC	10 ³ /μL	±2.2	21.2	21.2	20.54
RBC	10 ⁶ /μL	±0.25	5.99	5.99	6.08
HGB	g/dL	±0.7	19.5	19.5	19.57
HCT	%	±4.0	55.6	55.6	64.6
MCV	fL	±7.0	92.9	92.9	106.3
MCH	pg	±2.8	32.6	32.6	32.2
MCHC	g/dL	±3.0	35.1	35.1	30.3
PLT	10 ³ /μL	±70	531	531	549.3
MPV	fL	±3.0	5.1	7.9	8.3
RDW-CV *	%	±3.0	16.1	16.1	17.9
LY%	%	±10.0	21.2	21.2	14.83
MO%	%	±5.0	1.1	1.1	0.55
GR%	%	±10.0	77.7	77.7	84.62
LY	10 ³ /μL	±1.8	4.5	4.5	3.05
MO	10 ³ /μL	±0.7	0.2	0.2	0.11
GR	10 ³ /μL	±2.0	16.5	16.5	17.38

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.

Manual Registration

Lot	255
Level	High
Control	MEK-3D


<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU

TOKYO 161-8560 JAPAN



MEK-3DH



B255H



2025-08-05



0694-901686T

MEK-3DL

HEMATOLOGY CONTROL

2.0mL

ASSAY SHEET

LOT	B255L
	2025-08-05
OPEN	14DAYS

MEK-1301,1302,1305 Barcode Entry

The barcodes may be read in any order.

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
WBC	10 ³ /μL	±0.6	2.2	2.2	2.10
RBC	10 ⁶ /μL	±0.15	2.34	2.34	2.31
HGB	g/dL	±0.4	6.0	6.2	6.16
HCT	%	±2.0	18.6	18.6	20.7
MCV	fL	±9.0	79.3	79.3	89.4
MCH	pg	±2.4	25.6	26.5	26.7
MCHC	g/dL	±3.0	32.3	33.3	29.8
PLT	10 ³ /μL	±30	67	67	79.6
MPV	fL	±3.0	5.6	8.1	8.2
RDW-CV *	%	±3.0	18.2	18.2	16.6
LY%	%	±10.0	54.9	54.9	56.35
MO%	%	±5.0	7.6	7.6	4.06
GR%	%	±10.0	37.5	37.5	39.59
LY	10 ³ /μL	±0.5	1.2	1.2	1.18
MO	10 ³ /μL	±0.2	0.2	0.2	0.09
GR	10 ³ /μL	±0.5	0.8	0.8	0.83

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.

MEK-1305

QC ESR

45±17 mm

Manual Registration

Lot	255
Level	Low
Control	MEK-3D


<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU

TOKYO 161-8560 JAPAN



MEK-3DL



B255L



2025-08-05



0694-901684T