

MEK-3DN HEMATOLOGY CONTROL 2.0mL

ASSAY SHEET

LOT	B25YN
	2026-02-05
OPEN	14DAYS

MEK-1301,1302,1305 Barcode Entry

The barcodes may be read in any order. If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
WBC	10 ³ /μL	±0.8	7.6	7.6	7.66
RBC	10 ⁶ /μL	±0.20	4.82	4.82	4.76
HGB	g/dL	±0.5	13.9	14.1	14.12
HCT	%	±3.0	41.5	41.5	46.7
MCV	fL	±7.0	86.1	86.1	98.1
MCH	pg	±2.4	28.8	29.3	29.7
MCHC	g/dL	±3.0	33.5	34.0	30.2
PLT	10 ³ /μL	±50	240	240	257.1
MPV	fL	±2.0	5.3	8.0	8.0
RDW-CV *	%	±3.0	17.0	17.0	16.9
LY%	%	±10.0	35.1	35.1	32.82
MO%	%	±5.0	3.6	3.6	2.05
GR%	%	±10.0	61.3	61.3	65.13
LY	10 ³ /μL	±1.0	2.7	2.7	2.51
MO	10 ³ /μL	±0.4	0.3	0.3	0.16
GR	10 ³ /μL	±1.0	4.6	4.6	4.99

MEK-1305
QC ESR
9±5 mm

Manual Registration	
Lot	25Y
Level	Normal
Control	MEK-3D

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.



<https://www.nihonkohden.com/spi/index.html>



NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



MEK-3DN



B25YN



2026-02-05



0694-901685V

MEK-3DH	HEMATOLOGY CONTROL	2.0mL
----------------	---------------------------	--------------

ASSAY SHEET

LOT	B25YH
	2026-02-05
OPEN	14DAYS

MEK-1301,1302,1305 Barcode Entry

The barcodes may be read in any order.
 If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
WBC	10 ³ /μL	±2.2	21.3	21.3	20.89
RBC	10 ⁶ /μL	±0.25	5.98	5.98	5.95
HGB	g/dL	±0.7	19.2	19.2	19.37
HCT	%	±4.0	55.6	55.6	63.4
MCV	fL	±7.0	92.9	92.9	106.5
MCH	pg	±2.8	32.1	32.1	32.6
MCHC	g/dL	±3.0	34.5	34.5	30.6
PLT	10 ³ /μL	±70	515	515	537.6
MPV	fL	±3.0	5.4	8.1	8.4
RDW-CV *	%	±3.0	16.5	16.5	17.6
LY%	%	±10.0	14.5	14.5	12.43
MO%	%	±5.0	1.2	1.2	1.43
GR%	%	±10.0	84.3	84.3	86.14
LY	10 ³ /μL	±1.8	3.1	3.1	2.60
MO	10 ³ /μL	±0.7	0.3	0.3	0.30
GR	10 ³ /μL	±2.0	17.9	17.9	17.99

Manual Registration	
Lot	25Y
Level	High
Control	MEK-3D

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.



<https://www.nihonkohden.com/spi/index.html>


NIHON KOHDEN CORPORATION

 1-31-4 NISHIOCHIAI, SHINJUKU-KU
 TOKYO 161-8560 JAPAN


MEK-3DH



B25YH



2026-02-05



0694-901686V

MEK-3DL	HEMATOLOGY CONTROL	2.0mL
----------------	---------------------------	--------------

ASSAY SHEET

LOT	B25YL
	2026-02-05
OPEN	14DAYS

MEK-1301,1302,1305 Barcode Entry

The barcodes may be read in any order.
 If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
WBC	10 ³ /μL	±0.6	2.1	2.1	2.07
RBC	10 ⁶ /μL	±0.15	2.43	2.43	2.36
HGB	g/dL	±0.4	6.0	6.2	6.16
HCT	%	±2.0	18.7	18.7	20.8
MCV	fL	±9.0	77.0	77.0	88.0
MCH	pg	±2.4	24.7	25.5	26.1
MCHC	g/dL	±3.0	32.1	33.2	29.6
PLT	10 ³ /μL	±30	57	57	67.3
MPV	fL	±3.0	5.4	8.1	7.6
RDW-CV *	%	±3.0	17.5	17.5	15.7
LY%	%	±10.0	57.7	57.7	56.25
MO%	%	±5.0	10.0	10.0	6.67
GR%	%	±10.0	32.3	32.3	37.08
LY	10 ³ /μL	±0.5	1.2	1.2	1.16
MO	10 ³ /μL	±0.2	0.2	0.2	0.14
GR	10 ³ /μL	±0.5	0.7	0.7	0.77

MEK-1305
QC ESR
45±17 mm

Manual Registration	
Lot	25Y
Level	Low
Control	MEK-3D

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.



<https://www.nihonkohden.com/spi/index.html>


NIHON KOHDEN CORPORATION

 1-31-4 NISHIOCHIAI, SHINJUKU-KU
 TOKYO 161-8560 JAPAN


MEK-3DL



B25YL



2026-02-05



0694-901684V