


MK-RE2
HEMATOLOGY CONTROL
 for reticulocyte

3.0mL

ASSAY SHEET

LOT	MKRE265-2
	2026-07-05
OPEN	14 DAYS

		MEK-9200	RANGE
RET%	%	4.40	±2.60
RET	10⁶/μL	0.1549	±0.0915

2D Barcode Entry


• If RBC value is under 3.00 [$\times 10^6/\mu\text{L}$], please refrain from using.
 It may have deteriorated or hemolyzed, which affects the measurement accuracy.


<https://www.nihonkohden.com/spi/index.html>
Manual Registration

Lot	265					
Level	2					
Expiration	7	M	5	D	2026	Y
Control	MKRE					



 0123



NIHON KOHDEN CORPORATION
 1-31-4 NISHIOCHIAI, SHINJUKU-KU
 TOKYO 161-8560 JAPAN

MK-RE3	HEMATOLOGY CONTROL for reticulocyte	3.0mL
---------------	---	-------

ASSAY SHEET

LOT	MKRE265-3
	2026-07-05
OPEN	14 DAYS

		MEK-9200	RANGE
RET%	%	9.02	±2.98
RET	10⁶/μL	0.3181	±0.1051

2D Barcode Entry


• If RBC value is under 3.00 [$\times 10^6/\mu\text{L}$], please refrain from using.
It may have deteriorated or hemolyzed, which affects the measurement accuracy.



<https://www.nihonkohden.com/spi/index.html>

Manual Registration						
Lot	265					
Level	3					
Expiration	7	M	5	D	2026	Y
Control	MKRE					






	NIHON KOHDEN CORPORATION 1-31-4 NISHIOCHIAI, SHINJUKU-KU TOKYO 161-8560 JAPAN
---	--

MK-RE1

HEMATOLOGY CONTROL
for reticulocyte

3.0mL

ASSAY SHEET

 LOT	MKRE265-1
	2026-07-05
 OPEN	14 DAYS

		MEK-9200	RANGE
RET%	%	1.71	±1.29
RET	10 ⁶ /μL	0.0601	±0.0454

2D Barcode Entry



• If RBC value is under 3.00 [$\times 10^6/\mu\text{L}$], please refrain from using.
It may have deteriorated or hemolyzed, which affects the measurement accuracy.



<https://www.nihonkohden.com/spi/index.html>

Manual Registration

Lot	265					
Level	1					
Expiration	7	M	5	D	2026	Y
Control	MKRE					


0123





NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN