

MEK-3DN HEMATOLOGY CONTROL 2.0mL

ASSAY SHEET

LOT	B265N
	2026-08-05
OPEN	14DAYS

MEK-1301,1302,1305 Barcode Entry

The barcodes may be read in any order. If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
WBC	10 ³ /μL	±0.8	7.6	7.6	7.61
RBC	10 ⁶ /μL	±0.20	4.77	4.77	4.78
HGB	g/dL	±0.5	14.0	14.1	14.16
HCT	%	±3.0	41.4	41.4	47.1
MCV	fL	±9.0	86.7	86.7	98.6
MCH	pg	±2.4	29.4	29.6	29.6
MCHC	g/dL	±3.0	33.8	34.1	30.1
PLT	10 ³ /μL	±50	240	240	247.5
MPV	fL	±2.0	5.4	8.0	7.9
RDW-CV *	%	±3.0	16.8	16.8	17.1
LY%	%	±10.0	30.2	30.2	32.12
MO%	%	±5.0	3.1	3.1	3.21
GR%	%	±10.0	66.7	66.7	64.67
LY	10 ³ /μL	±1.0	2.3	2.3	2.44
MO	10 ³ /μL	±0.4	0.2	0.2	0.24
GR	10 ³ /μL	±1.0	5.1	5.1	4.93

MEK-1305
QC ESR
10±5 mm

Manual Registration	
Lot	265
Level	Normal
Control	MEK-3D

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.



<https://www.nihonkohden.com/spi/index.html>



NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



MEK-3DN



B265N




2026-08-05



0694-901685X

MEK-3DH	HEMATOLOGY CONTROL	2.0mL
----------------	---------------------------	--------------

ASSAY SHEET

LOT	B265H
	2026-08-05
OPEN	14DAYS

MEK-1301,1302,1305 Barcode Entry

The barcodes may be read in any order.
 If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
WBC	10 ³ /μL	±2.2	20.2	20.2	19.89
RBC	10 ⁶ /μL	±0.25	5.89	5.89	5.90
HGB	g/dL	±0.7	19.1	19.0	19.17
HCT	%	±4.0	55.0	55.0	62.8
MCV	fL	±9.0	93.3	93.3	106.5
MCH	pg	±2.8	32.4	32.3	32.5
MCHC	g/dL	±3.0	34.7	34.5	30.5
PLT	10 ³ /μL	±70	499	499	497.1
MPV	fL	±3.0	5.3	7.9	8.1
RDW-CV *	%	±3.0	16.3	16.3	17.8
LY%	%	±10.0	15.8	15.8	13.30
MO%	%	±5.0	1.0	1.0	1.15
GR%	%	±10.0	83.2	83.2	85.55
LY	10 ³ /μL	±1.8	3.2	3.2	2.65
MO	10 ³ /μL	±0.7	0.2	0.2	0.23
GR	10 ³ /μL	±2.0	16.8	16.8	17.01

Manual Registration	
Lot	265
Level	High
Control	MEK-3D

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.



<https://www.nihonkohden.com/spi/index.html>



NIHON KOHDEN CORPORATION

1-31-4 NISHIOCHIAI, SHINJUKU-KU
 TOKYO 161-8560 JAPAN



MEK-3DH



B265H




2026-08-05



0694-901686X

MEK-3DL HEMATOLOGY CONTROL 2.0mL

ASSAY SHEET

LOT	B265L
	2026-08-05
OPEN	14DAYS

MEK-1301,1302,1305 Barcode Entry

The barcodes may be read in any order. If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
WBC	10 ³ /μL	±0.6	1.8	1.8	1.81
RBC	10 ⁶ /μL	±0.15	2.42	2.42	2.36
HGB	g/dL	±0.4	5.9	6.0	6.04
HCT	%	±2.0	18.2	18.2	19.7
MCV	fL	±9.0	75.2	75.2	83.6
MCH	pg	±2.4	24.4	24.8	25.6
MCHC	g/dL	±3.0	32.4	33.0	30.7
PLT	10 ³ /μL	±30	58	58	65.9
MPV	fL	±3.0	6.0	8.3	8.2
RDW-CV *	%	±3.0	18.1	18.1	16.8
LY%	%	±10.0	54.3	54.3	59.13
MO%	%	±5.0	6.5	6.5	3.21
GR%	%	±10.0	39.2	39.2	37.66
LY	10 ³ /μL	±0.5	1.0	1.0	1.07
MO	10 ³ /μL	±0.2	0.1	0.1	0.06
GR	10 ³ /μL	±0.5	0.7	0.7	0.68

MEK-1305
QC ESR
48±17 mm

Manual Registration	
Lot	265
Level	Low
Control	MEK-3D

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



MEK-3DL



B265L



2026-08-05



0694-901684X